



## Security by Password Protection

Secura® comes with password protection to restrict adjustment of the balance settings only to authorized persons. If you set a password, all functions that could change the metrological performance specifications are inactivated.

## Interactive User Guidance for Fast GLP- | GMP-Compliant Printouts

If you require traceable documents, Secura® will guide you through the steps for obtaining GLP- | GMP-compliant printouts. The novel user guidance will support you especially for working with traceable sample identification.

Used in conjunction with the new premium YDP30 laboratory printer, which is automatically detected by Secura®, the balance generates printouts of the results on standard printer paper or on continuous labels.

2013-02-21	14:46
Sartorius	
Mod.	SECURA224-1S
SerNo.	0027703716
BAC:	00-50-02
APC:	01-70-02
Comp1 +	1.7082 g
Comp2 +	1.5017 g
Comp3 +	0.0510 g
T-Comp+	3.2608 g
<hr/>	
2013-02-21	14:47
Name:	